



QuantumDye™ Terminator v3.1 Cycle Sequencing Kit

DESCRIPTION

The QuantumDye™ v3.1 Terminator Cycle Sequencing Kit is based on the Sanger Dideoxynucleotide Chain Termination method. The kit is provided with a 2.5x concentrated ready-reaction premix, plus sequencing dilution buffer and controls (plasmid DNA template and sequencing primer) fully optimized for a highly flexible chemistry, designed for all kinds of sequencing applications, including *de novo* sequencing and resequencing. The kit generates data with uniform peak heights and optimized signal balance to produce long, high-quality reads.

KIT CONTENTS


QSD3-1000 QuantumDye Terminator Cycle Sequencing Kit v3.1; 1,000 reactions:

	Component	Concentration	Quantity
•	QuantumDye Cycle Sequencing Mix	2.5x	10 x 800 µl
•	QuantumDye Sequencing Buffer	5x	1 x 12 mL
•	pGEM® Control DNA	100 ng/µl	1 x 250 µl
•	M13(-21) Control Primer	5 µM	1 x 200 µl

PROTOCOL

The QuantumDye Terminator Kit was developed as a direct, drop-in substitute for BigDye® Terminator v3.1, with no changes in protocol, dilution, workflow, calibration or settings using DyeSet Z. Simply dilute the QuantumDye Cycle Sequencing Mix and substitute it for BigDye Terminator as described in your SOP. Using 0.5 - 1.0 µl of QuantumDye Cycle Sequencing Mix (1:16 - 1:8 dilution) produces results that are identical to BigDye Terminator v3.1 in comparison reactions with uniform peak heights, extremely low background and long sequence reads with high confidence scores.

©2020 QuantumSeq. All rights reserved. QuantumSeq, QuantumSeq logo and QuantumDye are trademarks of QuantumSeq. pGEM is a registered trademark of Promega Corporation. BigDye is a registered trademark of Thermo Fisher Scientific or its subsidiaries in the U.S. and certain other countries.



QuantumSeq

For research use only. Not for use in diagnostic procedures.

Made in USA